



Side Wedges with Connecting Sheet Fitting Instructions



Ensure the bed is at a safe and comfortable working height for yourself.



1. Place the connecting sheet on the bed, blue side up.



2. Pass both straps under the mattress and connect the buckles up on the other side. Adjust to tighten.



3. Secure a wedge on one side by attaching the hook and look strips, ensuring the right angle of the wedge is aligned to the edge of the bed. Repeat on the other side.

When using any bed accessories, make sure you have done your own risk assessment or follow the guidance of someone who has done a risk assessment. Check all equipment for signs of wear and tear prior to use.

2022-v. 1



MIP

The Softer Side of Healthcare