



 Features & Benefits

Our high density foam filled fall-out mats are designed for use with low level beds to ensure that patients do not injure themselves if they were to fall out of a bed with no rails. Made with a blue pvc, wipe clean material the fall-out mat is very easy to store and maintain and comes complete with durable carry straps. Flame Retardant to British Standard BS 7177.

Cleaning: Wipe Clean only. Wipe down with a mild detergent or a sterilising fluid. Do not use bleach or biological cleaners.

 Specifications

STYLE	DESCRIPTION	FABRIC CONTENT	SIZE
139/W/PU	Fall-Out Mats	PVC with foam filling	60 x 190 x 5cm