



Features & Benefits

Our high density foam filled pair of side wedges with a breathable connecting sheet are designed to prevent residents from falling out of beds with no rails. The side wedges are attached to the connecting sheet using a large hook and loop strip, which is then attached underneath the bed using durable straps and quick release clips. Flame Retardant to British Standard BS 7177.

Cleaning: Wipe Clean only. Wipe down with a mild detergent or a sterilising fluid. Do not use bleach or biological cleaners.

Specifications

| STYLE | COLOUR/PATTERN | DESCRIPTION | FABRIC CONTENT | SIZE |
|------------|----------------|-----------------------------------|---|-----------|
| 138/W/PU | Blue | Side Wedges with Connecting Sheet | Foam filled PVC wedges with PU connecting sheet | 20 x 90cm |
| 138/WEDGES | Blue | Side Wedges | Foam filled PVC wedges | 20 x 90cm |
| 137/W/PU | Blue | Connecting Sheet | PU connecting sheet | 90 x 90cm |