



MIP (UK) Ltd
2nd Floor, Clipper House
Leighton Industrial Park
Billington Road
Leighton Buzzard
Bedfordshire, LU7 4AJ

Tel: 01525 218146
Fax: 01525 218147

Email: mipuk@mipinc.com
Web: www.mipuk.co.uk

Patient Hoist Sling **Compatibility Statement**

It is now widely acknowledged that slings can be interchanged with hoists provided both hoist & sling are fully and correctly assessed, inspected and maintained. The spreader bar must be compatible with the sling that is to be used, i.e. loop sling to straight/coat hanger type spreader bar or keyhole/clip type to stud fixing spreader bar.

Cromptons Healthcare Slings supplied by MIP UK Limited have been tested and meet the requirements on the following hoists using either a clipped cradle, a standard loop fitting coat hanger or 'T' Spreader bar arrangement:

Oxford, Arjo, Liko, Invacare, Guldmann, Steiss, Oxford Presence, Parker Bath Hoist, Wispa Duo, Porta 250, Liko Viking M, Liko Glovo, Liko Uno 102, Liko Sabina li, Arjo Walker, Arjo Encore, Arjo Maxi Move, Arjo Maxi Sky, Arjo Sara Plus, Arjo Sara 2000, Molift Partner

In the event that the hoist you are working with does not appear on the list, please contact us and we will be more than happy to carry out a compatibility assessment.

Public Liability Insurance

Our insurers provide Public and Products liability cover to the sum of ten million pounds in respect of injury or damage in the use of slings supplied by our company. Full insurance certificates are available upon request.

Contact Information

MIP (UK) Limited
2nd Floor, Clipper House
Leighton Industrial Park
Billington Road
Leighton Buzzard
Bedfordshire
LU7 4AJ

Tel: 01525 218146
Fax: 01525 218147

Email: mipuk@mkipinc.com
Web: www.mipuk.co.uk



Registered in England: 2695429