



FUSION™ Overview

Discover the Fusion® Washable Bed Pad by MIP – the latest in bed pad technology.

The Fusion® Bed Pad is designed with both comfort and protection in mind, offering advanced skin care for those dealing with incontinence. Thoughtfully crafted with a fused, un-quilted design, Fusion® helps maintain skin integrity by wicking moisture away from the body. This keeps the care recipient cool, dry, and comfortable whilst lying on the pad.

Fusion® features a smooth surface and a specially engineered brushed bottom barrier ensuring the pad stays in place. This minimises the risk of the pad creasing and rucking on the bed, therefore avoiding pressure care issues. By combining innovation with compassionate care, Fusion® supports both caregivers and care receivers, ensuring a restful and worry-free experience.

MIP Summary

As a global supplier and engineering expert of textile solutions, MIP provides comfort and efficiency for the healthcare community and its people.

We believe in quality without compromise. Even the smallest improvements can hold the greatest importance.

At MIP, we have a wide textiles product offering including ranges in Bedding, Linen, Laundry, Bags, Linen Collection & Distribution, Apparel, Continence Care and Moving & Handling.


MIP. The softer side of healthcare.



For more information, visit our website at www.mipuk.co.uk

Specifications:

Dimensions	91 × 91cm
Colour	Turquoise/Blue Barrier
Absorbency	2 Litre
Description	'This Side Down' Print Latex Free Materials



Features:



Moisture Wicking
Helps the care receiver's body remain cool and increasing comfort



Fused Construction
Fused layers makes the bed pad more resilient and efficient to launder



Superior Absorbency
The soaker layer has a 2L absorbency



Brushed Barrier
Provides protection to the bedding and mattress, as well as keeping the pad in place



Reduce Waste
Can be washed multiple times unlike disposables, reducing landfill waste



Supports Pressure Care
Smooth surface minimises the risk of creases to avoid pressure care issues



Less Bed Changes
As Fusion® stays in place, there is less risk of leaks on the sheets, meaning less bed changes are needed



Thermal Disinfection
Fusion® can be washed at thermal disinfection temperatures

Washing Instructions:

 Wash **up to 95 degrees**
 Tumble on **medium heat**

 **Do not iron**
 **Do not bleach**

