



Patient Specific Tubular Flat



## Features & Benefits

**MIP Slide Sheets** are made from Transtex, a specially developed fabric that is silicone coated on both sides to provide a slick, low friction finish for easy patient movement and turning.

Two-Tone Tubular Slide Sheets are formed by sewing the fabric into a continuous tube with a double sewn low profile seam that reduces potential tissue viability issues whilst adding strength and durability. The unique two-tone design helps simplify correct product orientation and can also support training and best practice during manual handling procedures.

A web hanging loop is provided to allow the tube to be kept with the patient at all times.

For all repositioning and manual handling solutions it is recommended to seek advice / approval from a manual handling expert.

CE and CA Marked

-----

Available to purchase through NHS Supply Chain [here](#).

Explore our Slide Sheet Blogs for expert guidance and practical insights:

- ? [Slide Sheet FAQs](#)
- ? [Flat vs Tubular Slide Sheets](#)
- ? [Why Choose a Reusable Slide Sheet?](#)
- ? [Caring for your Slide Sheets](#)
- ? [Checking your Slide Sheets](#)
- ? [The Benefits of using Slide Sheets](#)

## Specifications

| STYLE        | COLOUR/PATTERN | DESCRIPTION                                         | SIZE        |
|--------------|----------------|-----------------------------------------------------|-------------|
| PSTL200140AL | Aqua and Lilac | Tubular Patient Use Bariatric Two-Tone Slide Sheets | 140 x 200cm |